

## SPECIFICATIONS :

### Characteristics :

**Mode of operation :**  
Bi-level, constant pressure

**Pressure output :**  
**EPAP :** from 3 cm H<sub>2</sub>O to IPAP or 20 cm H<sub>2</sub>O (+/-0.5 cm H<sub>2</sub>O)  
**IPAP :** from EPAP or 3 cm H<sub>2</sub>O to 25 cm H<sub>2</sub>O (+/-0.5 cm H<sub>2</sub>O)  
**Comfort pressure :** from 3 cm H<sub>2</sub>O to EPAP

**Ramp duration :** 0 to 45 minutes adjustable by 5 minutes

**Advanced setting parameters :**  
ISENS & ESENS : from 1 to 10  
Pressure Rise Time : from 1 to 5  
Back-up rate (DuoST) : from 0 to 25 bpm  
I/E ratio (DuoST) : from ¼ to 1/0.5

**Max. flow for 1 cmH<sub>2</sub>O pressure drop at the patient output :**  
180 lpm at 20 cmH<sub>2</sub>O without humidifier

**Noise level (with reservoir):** 27 dB measured with an uncertainty of 2 dB according to ISO 17510-1:2002

**Humidifier chamber capacity :**  
300 ml at the maximum water level

**Heater adjustable range and wattage :**  
0 to 10 adjustable (40 Watts at setting 10)

### Physical Characteristics :

**Size :** 20 cm wide X 23 cm deep X 12 cm high

**Weight :** 1,5 kg with humidifier; 1,4 kg without humidifier

**Filter :** Cleanable standard or disposable HEPA

### Power characteristics :

**Power supply :** Input range 100-240 V, 50-60 Hz

**Power supply on battery :** 13 V

**Power consumption :**  
50 VA max without humidifier; 90 VA with humidifier at setting 10

**IEC 60601-1 classification :**  
Class II (double isolation) Type BF

### Environmental requirements :

**Operating temperature**  
5°C to 35°C, 10% to 95% relative humidity (RH), non-condensing

**Storage temperature :**  
-20°C to 60°C, relative humidity (RH) ≤95%, non-condensing

**Operating altitude :** from sea level to 2700 meters

## ORDERING INFORMATIONS

### Products and accessories :

DreamStar™ Duo (includes patient circuit 1,8 m and carrying case)	M-115503- UK
DreamStar™ Duo Evolve w/reservoir (includes patient circuit 1,8 m and carrying case)	M-115523- UK
DreamStar™ Duo Evolve w/ cover (includes patient circuit 1,8 m and carrying case)	M-115513- UK
DreamStar™ DuoST (includes patient circuit 1,8 m and carrying case)	M-115504- UK
DreamStar™ DuoST Evolve w/reservoir (includes patient circuit 1,8 m and carrying case)	M-115524- UK
DreamStar™ DuoST Evolve w/ cover (includes patient circuit 1,8 m and carrying case)	M-115514- UK
DreamStar™ Water Reservoir (includes base, cover and gasket)	M-315510-14
DreamStar™ Cover	M-315530-19
Filter holder kit (27x42)	M-314820-02
96 Inlet Air filter (27x42)	M-414840-06
100 High efficiency air filter (27x42) Carrying case DreamStar™	M-414841-07 M-815505-00
Power cord 2 poles - UK	M-660601-40
Cigarette lighter cable 3A	M-214530-01
Battery cable 3A ring term	M-214530-00
Patient circuit 1,80 m / diam 15mm	M-261000-00
Patient circuit 1,80 m / diam 22mm	M-261000-04
Patient manual DreamStar™ Duo / DuoST	M-155DFU00-30
Practitioner manual DreamStar™ Duo / DuoST Evolve	M-155DFU00-70

### Data Management :

DreamStar™ Analyze software	M-215630-00
DreamStar™ Analyze clinical kit	M-215630-01
Memory card DreamStar™	M-215530-00
Memory card reader	M-214830-03
User manual	M-156DFU00-20

### References :

M-115503- UK

M-115523- UK

M-115513- UK

M-115504- UK

M-115524- UK

M-115514- UK

M-315510-14

M-315530-19

M-314820-02

M-414840-06

M-414841-07

M-815505-00

M-660601-40

M-214530-01

M-214530-00

M-261000-00

M-261000-04

M-155DFU00-30

M-155DFU00-70

M-215630-00

M-215630-01

M-215530-00

M-214830-03

M-156DFU00-20




# DreamStar™ Duo & Duo Evolve DreamStar™ DuoST & DuoST Evolve

## Bi-Level Devices



DreamStar™ Duo


 The very latest innovations of the DreamStar™ Duo & DreamStar™ DuoST support optimum function of these devices used for non-invasive ventilation of patients suffering from Sleep Apnea Syndrome or respiratory failure but who are not dependant on ventilatory assistance.

These technological innovations are designed to maximize comfort and patient compliance with the therapy allowing an easy treatment adaptation and follow up.

Both DreamStar™ Duo & DuoST detect residual apneas and hypopneas, provide the estimated tidal volume and the leak level to easily evaluate the effectiveness of

## CONFIGURATIONS

The DreamStar™ Duo & DreamStar™ DuoST have a configuration without humidifier and the DreamStar™ Duo & DreamStar™ DuoST Evolve with an integrated heated humidifier. The Evolve configurations may be used either with a water reservoir for humidification or without humidification by replacing the water reservoir with a cover.

 Quick and simple setting of the humidifier from the main screen of the device.




DreamStar™ DuoST Evolve w/ water reservoir

For more information contact us at the following addresses :  
customerservice@sefam-medical.com  
marketing@sefam-medical.com


DreamStar™ is a trademark of SEFAM

CE 0459



 **The AUTO ON function** enables patients to begin their treatment by simply breathing into their mask. They still have the option of activating the device by pressing the on/off button.

**Specific AUTO-CALIBRATION function** ensures that the prescribed pressure is delivered at the patient's interface. This makes the device compatible with almost any combination of mask and patient circuit of 22mm or 15mm diameter.

 **Settings ADAPT feature.** The activation of this function, when the device is operational, allows to display a page of main screen indicating the values of tidal volume, respiratory rate, average leak and inspiratory time of the last cycle. This feature allows the practitioner to define optimum treatment parameters in a simpler and fastest way.



**Flow Sens Algorithm** The unique "FlowSens" feature allows Dreamstar™ Duo & Duo ST devices to adjust automatically the triggering levels of inspiratory and expiratory phases of the patient breathing and to answer instantaneously to any patient respiratory changes. FlowSens detects and compensates mask or patient circuit leaks but also adjust the device parameters to ensure optimal therapy performances in presence of leaks. In addition, to be able to answer to other special patient needs, the advance setting menu gives access to customized setting for these inspiratory and expiratory trigger levels and the pressure rise time settings.

#### COMPLIANCE DATA Management

Compliance data are stored for 365 days in the internal memory of the device. The last 16 hours of the patient's actual flow data provides in addition extra useful information about the effectiveness of the therapy.

The **DreamStar™ Analyze software** is a user friendly software which enables downloading, analyzing and printing the compliance data to a fully configurable report. The data can be retrieved either entirely by using the mini USB DreamStar™ cable or through the DreamStar™ memory card for the last 40 eight-hour sessions.

#### Dreamstar™ water & humid control technology

The innovative design of the water reservoir "DreamStar™ Water Control" and the automatic humidification system "Humid Control" provide a new concept in water management. This prevents the risk of water return to the devices and minimizes condensation in the patient circuit to maintain a stable level of humidification overnight, for a therapy combining comfort and optimal efficiency.

#### FUNCTIONALITIES :

##### Ease of use :

- Front-facing screen
- Auto reminder settings
- Large, backlit buttons
- Multilingual menu
- Can run on external battery (by-pass humidification option)
- Download compliance data via DreamStar™ card or mini USB cable
- Upload settings via the DreamStar™ card or mini USB cable
- Complete software analysis DreamStar™ Analyze

##### Technologies :

- Auto ON function
- Settings ADAPT feature
- Specific Auto-calibration function
- FlowSens Algorithm
- DreamStar™ Water & Humid Control+ technology in the Evolve configuration

##### Fully integrated :

- Integrated HC humidifier
- Single electrical connection for both device and heated humidifier

##### Automated :

- Automatic altitude and barometric pressure change compensation
- Universal internal power supply
- Automatic leak compensation

#### COMPACT DESIGN

DreamStar™ Duo & DuoST in the Evolve configuration have a fully integrated heated humidifier for enhanced breathing comfort, while remaining one of the most compact devices in the market. The humidifier is located at the rear of the devices with its water reservoir. The reservoir and the air inlet filter are easy to remove for cleaning.

The fully integrated power supply and heated humidifier makes DreamStar™ Duo & DuoST Evolve one of the lightest devices in the market. This ensures maximum comfort while using the device, either when travelling or being at home.