



Actual Size: 140 mm (width)
x 175 mm (height)

Ruler is for scaling purposes.

How to use the fitting template:


1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Follow directions on template to determine appropriate cushion size

Mirage™ FX

Mirage™ FX for Her

Mirage is a trademark and/or registered trademark of the ResMed family of companies.
© 2015 ResMed Ltd. 1018996/1 2015-08

Rx Only  0123



Sample size check

Mirage™ FX


Mirage™ FX for Her

Small

Standard

Wide

ResMed

Rx Only  0123

© 2015 ResMed Ltd. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.